

Model : Nitroplus - 15 C



ADVANTAGES

- 15 place Block COD Digestion system
- Touch screen operation
- Space saving, compact and user friendly.

MAIN FEATURES

- Fifteen place COD digestion system
- Casted Aluminium alloy heater
- High grade Extruded aluminium block
- Insulation by high grade ceramic wool
- Acid reflux by Glass air condenser
- Stainless steel rack for handling digestion tubes
- Mild steel support rack to place air condensers
- Plastic tray with a provision to add reagents to the sample tube in water
- Flat bottom digestion tube enabling titration of sample without transfer

TECHNICAL DATA

Description	Nitroplus - 15CE	Nitroplus - 15CT
No of samples	15	15
Tube dimension	36 X 200 mm	36 X 200 mm
Tube capacity	150 ml	150 ml
Process controller	7 segment LED controller	4.3" TFT touch panel
Heater wattage	2500w	2500w
Temperature / Time program	-	5 program / 5 steps each
Temperature range	50° C to 300 °C	50° C to 300° C
Temperature accuracy	+/- 1° C	+/- 1° C
Temperature sensor	K type thermocouple	K type thermocouple
Insert rack	Stainless steel	Stainless steel
Electrical requirement	230VAC / 50Hz / Single phase	230VAC / 50Hz / Single phase

Manufactured by :

Pelilab Scientific Private Limited

No. 4, Anna Street, Srinivasa Nagar, Padi, Chennai - 600 050.
 Web : www.pelilabscientific.com Email : info@pelilabscientific.com

Marketed by :