

Model : D waterplus



ADVANTAGES

- Automatic water distillation system
- Easy dismantling of water condenser for easy cleaning.
- Optical sensor for maintaining the water level.
- Space saving, compact and user friendly.
- Double Distillation unit by integrating with another same unit.

MAIN FEATURES

- Stainless steel steam generator
- Auto refill of water during distillation.
- Condenser water saving mechanism
- Water to the steam condenser is cut off when sufficient water is collected in the reservoir.
- Can be integrated to second distillation system to collect double distilled water
- Collection of distilled water in 5 lit container
- Automatic switching ON of system when water level in the container is below the fixed level.
- Membrane keypad with LED indications for Power, Heater and water condenser.

TECHNICAL DATA

Description	D waterplus
Steam generator	High grade stainless steel
Water condensor	High grade stainless steel
Water Level monitoring and control	By Optical sensor
Water reservoir	PVC
Reservoir capacity	5 lit
Rate of distillation	3.5 lit/Hr.
Heater	Immersion heater
Heater wattage	2000 w
Electrical requirement	230 v/50Hz/Single Phase

Manufactured by :

Pelilab Scientific Private Limited

No. 4, Anna Street, Srinivasa Nagar, Padi, Chennai - 600 050.
 Web : www.pelilabscientific.com Email : info@pelilabscientific.com

Marketed by :