

## Model : Wamerplus



### ADVANTAGES

- Hotair oven and Incubator in single unit.
- Touch screen operation
- Time and Temperature programmable operation.
- Optional PC connectivity and Data logging.

### MAIN FEATURES

- Forced air circulation by convection
- Auto cut off of heater & when door is opened.
- Perfect insulation with high grade ceramic wool.
- 4.3" TFT touch operation
- Feather touch membrane key pad
- Over temperature protection by thermostat.
- Inner stainless steel chamber, outer mild steel with powder coated finish.
- Round folding at the four corners for easy cleaning.
- Tray height adjustable to 7 heights

### TECHNICAL DATA

Description	HCI-91E	HCI-91T	HCI-91TP
Capacity	91 ltr	91 ltr	91 ltr
Inner dimension	450 x 450 x 450 mm	450 x 450 x 450 mm	450 x 450 x 450 mm
Outer dimension	690 x 530 x 530 mm	690 x 530 x 530 mm	690 x 530 x 530 mm
Process controller	7 segment LED controller	4.3" TFT touch panel	4.3" TFT touch panel with PC connectivity
Heater wattage	Oven - 1500w, incubator - 500w	Oven - 1500w, incubator - 500w	Oven - 1500w, incubator - 500w
Temperature / Time program	No	5 Program/5 steps each	5 Program/5 steps each
Temperature range	5° above ambient to 200° C	5° above ambient to 200° C	5° above ambient to 200° C
Temperature accuracy	+/- 0.1° C	+/- 0.1° C	+/- 0.1° C
Temperature uniformity	+/- 2° C	+/- 2° C	+/- 2° C
PC connectivity & Data Logging	No	No	Yes
No of racks	2	2	2
Electrical requirement	230 v/50Hz/Single Phase	230 v/50Hz/Single Phase	230 v/50Hz/Single Phase

Manufactured by :

**Pelilab Scientific Private Limited**

No. 4, Anna Street, Srinivasa Nagar, Padi, Chennai - 600 050.  
 Web : [www.pelilabscientific.com](http://www.pelilabscientific.com) Email : [info@pelilabscientific.com](mailto:info@pelilabscientific.com)

Marketed by :