



Pelilab Scientific Private Limited

AUTOMATIC WATER DISTILLATION SYSTEM

D-WATERPLUS-1E

The D-Waterplus-1E is a compact automatic water distillation system designed for labs requiring a steady supply of pure water

With a stainless steel steam generator & condenser, 5L storage tank, and auto shut-off controls, it delivers safe, reliable, and maintenance-friendly operation.

Key Benefits:

- ✓ Automatic operation water-level control
- ✓ Compact & space-saving design
- ✓ Auto cut-off & overload protection
- ✓ Easy to clean & maintain
- ✓ Expandable to double distillation



TECHNICAL SPECIFICATIONS :

- Distillation Rate: 2.5–3 L/hr
- Construction: Stainless steel steam generator & condenser
- Heater: 2500 W immersion heater
- Storage Tank: 5L capacity
- Controls: Membrane keypad with indicators
- Safety: Auto cut-off, over-voltage protection, water-level monitoring
- Operation: Fully automatic ON/OFF by water level
- Cleaning: Easy dismantling for scale removal
- Expandable: Convertible to double distillation system
- Power: 230V / 50Hz / 16A

APPLICATIONS:

- Clinical & Research Laboratories
- Pharmaceutical & Chemical Labs
- Microbiology & Biotechnology Labs
- Quality Testing Facilities



Pelilab Scientific Private Limited



Visit the Website

PELILAB SCIENTIFIC PRIVATE LIMITED

☎ +91 99655 17164

✉ info@pelilabscientific.com

🌐 www.PeliLabScientific.com

📍 No.4, Anna Street, Srinivasa Nagar,

📍 Akshaya Colony, Padi,

Chennai - 600 500,

Tamil Nadu, India.