



Pelilab Scientific Private Limited

REFRIGERATED INCUBATOR SHAKER

ORBIPLUS-ET

The ORBIPLUS-ET combines orbital shaking, precise heating, and integrated refrigeration in a compact, high-performance unit.

Designed for microbiology, biotechnology, and cell culture applications, it ensures accurate temperature control (5°C–65°C), noiseless operation, and user-friendly programming through a 4.3" TFT touchscreen controller.

Applications:

- ✓ Microbiology & Biotechnology Research
- ✓ Cell Culture & Fermentation Studies
- ✓ Pharmaceutical & Food Industry R&D
- ✓ Environmental & Agricultural Testing Labs
- ✓ Universities & Research Institutes



TECHNICAL SPECIFICATIONS :

- Orbital motion: 50–400 RPM, $\pm 1\%$ accuracy
- Temperature range: 15°C to 65°C ($\pm 0.1^\circ\text{C}$)
- 4.3" TFT touchscreen with auto & manual modes
- Programmable: 5 programs $\times$ 5 steps (RPM, Temp, Time)
- In-built compressor cooling
- Dual removable anodized aluminium trays (460 $\times$ 460 mm)
- Universal tray fits 250 ml–2000 ml flasks, 25 clamps included
- Insulated hood & acrylic viewing door with auto safety cut-off

KEY BENEFITS:

- ✓ ISO Certified Company
- ✓ Wide Temperature Range (5°C – 65°C)
- ✓ Built-in Refrigeration System
- ✓ Programmable Multi-Step Operation
- ✓ Energy-Efficient & Noiseless Design
- ✓ Trusted Pelilab After-Sales Support
- ✓ Integrated fluorescent & UV illumination



Pelilab Scientific Private Limited



Visit the Website

PELILAB SCIENTIFIC PRIVATE LIMITED

☎ +91 99655 17164

✉ info@pelilabscientific.com

🌐 www.PeliLabScientific.com

📍 No.4, Anna Street, Srinivasa Nagar,

📍 Akshaya Colony, Padi,

Chennai - 600 500,

Tamil Nadu, India.