



Pelilab Scientific Private Limited

HOT AIR OVEN CUM INCUBATOR

WAMERPLUS-091T

The WAMERPLUS-091T is a versatile laboratory solution that combines incubator and the hot air oven functionalities in a single unit.

Designed with precision temperature control, advanced safety features, and user-friendly operation, it ensures reliable performance for microbiology, pharmaceutical, research, and quality testing applications.

Applications:

- ✓ Microbiology Labs
- ✓ Pharmaceutical R&D
- ✓ Industrial QC Testing
- ✓ Biotechnology & Food Safety Research



KEY HIGHLIGHTS:

- Dual Operation – Hot Air Oven & Incubator in one unit
- Precise PID Temperature Control (5°C above ambient to 150°C)
- ISO-grade Safety – Over-temperature & Over-voltage protection
- 4.3” TFT Touchscreen Controller – Smart, user-friendly interface
- Forced Air Circulation – Uniform heating & incubation
- Robust MS outer body with stainless steel inner chamber
- Backed by Pelilab’s trusted after-sales service

DETAILED SPECS:

- Capacity: 91 Liters
- Inner Chamber: 450 x 450 x 450 mm
- Temperature Range (Incubator): 5°C above ambient to 70°C
- Temperature Range (Oven): 5°C above ambient to 150°C
- Accuracy: ±0.5°C at 37°C, ±1% at 120°C
- Heater Wattage: Oven – 1625W | Incubator – 325W
- Trays: 2 adjustable (7 height levels)
- Safety: Auto cut-off (door open), Thermostat protection, Circuit breaker
- Other Features: Glass viewing window, Easy-to-clean rounded corners, Additional port for temp cross-check



Pelilab Scientific Private Limited



Visit the Website

PELILAB SCIENTIFIC PRIVATE LIMITED

+91 99655 17164
info@pelilabscientific.com
www.PeliLabScientific.com
No.4, Anna Street, Srinivasa Nagar,
Akshaya Colony, Padi,
Chennai - 600 500,
Tamil Nadu, India.